

Beeps, Buzzers and Alarms

The Hemodialysis machine

Madhukar Misra MD
Professor of Medicine

University of Missouri Columbia

The systems Approach

- The HD machine cannot be isolated from the rest
- It is a part of an assembly
 - **The HD machine**
 - The extracorporeal circuit
 - The dialyzer
 - The water treatment system
 - Dialysate concentrate supply
 - The data network

Is there an Ideal Hemodialysis Machine?

- The ‘ideal’ machine does not exist and will never exist
- The utilitarian philosopher’s ‘ideal’ :
 - Achieve minimum dialysis dose
 - Achieve accuracy in controlled values
 - Achieve quality while controlling costs
- Machines are smarter in sensing/measuring/balancing
- User errors are common.

**SEARCH**

Most Popular Searches

[Home](#) [Food](#) [Drugs](#) [Medical Devices](#) [Radiation-Emitting Products](#) [Vaccines, Blood & Biologics](#) [Animal & Veterinary](#) [Cosmetics](#) [Tobacco Products](#)

MAUDE - Manufacturer and User Facility Device Experience

[FDA Home](#) [Medical Devices](#) [Databases](#)

JAMA

The Journal of the
American Medical Association

[Home](#) [Current Issue](#) [All Issues](#) [Online First](#) [Collections](#) [CME](#) [Multimedia](#)

June 2, 1993, Vol 269, No. 21 >

< Previous Article

Full content is available to subscribers
[Subscribe/Learn More](#)

Next Article >

ARTICLE | June 2, 1993

Introducing MEDWatch A New Approach to Reporting Medication and Device Adverse Effects and Product Problems

David A. Kessler, MD; Sharon Natanblut, MPA; Dianne Kennedy, MPH, RPh; Eliot Lazar, MD; Peter Rheinstein, MD, JD, MS; Chuck Anello, ScD; Dave Barash, RPh; Ilisa Bernstein, PharmD; Ross Bolger, RPh; Kay Cook, JD; Mary Pat Couig, RN, MPH; Jerry Donlon, MD, PhD; Joyce Johnson, DO, MA; Catherine Lorraine, JD; Tom McGinnis, RPh; John Nazario, RPh; Stuart Nightingale, MD; Carl Peck, MD; Mary Pendergast, JD; Suresh Rastogi, PhD; Chet Reynolds, MBA; Renie Schapiro, MPH; Linda Tollefson, DVM; Ann Wion, JD

JAMA. 1993;269(21):2765-2768. doi:10.1001/jama.1993.03500210065033.

Text Size: A A A

Machine malfunction is rare

Hemodialysis is prone to user errors

- Poor fixation of cannulas
- Incorrect setting of alarm limits and control values

Failure to follow Protocol is common

Event type	Number	% of total events, N=526
Medication error	150	28.5
Failure to follow protocol	68	12.9
Lab/Blood bank related	52	9.9
Procedure complication	45	8.6
Needle disconnection	32	6.1
Needle infiltration	32	6.1
Equipment failure	25	4.8
System Clotting	23	4.4


Fundamental Requirements of Hemodialysis Instruments

Practical Aspects

- Propel blood & dialysate
- Provide life sustaining “dose”
- Anticoagulation
- Maintain body temperature
- Prime, rinseback, and rapid hydration
- Add meds/draw blood

Hazard Prevention

- Blood loss to environment
- Hemolysis
- Air embolism
- Excessive or inadequate fluid removal
- Infection
- Acute or chronic toxicity
- Anaphylaxis



Beeps, Buzzers and Alarms

Blood Circuit



Inflow pressure

Outflow pressure

Air Detector

Dialysate Circuit




Conductivity

Temperature

Hemoglobin

Safety monitors

- Can not ‘prevent’ adverse events
- Can only ‘detect and mitigate’ the potential harm



Adapted from Handbook of Dialysis Therapy. Nissenson and Fine. 4th edition,
Saunders Elsevier Pubs.


Know your machine!

No Alarm is Fail-Safe

- Alarms represent potential malfunction in the system
- NEVER alter alarm sensitivity and range
- Discourage Alarm disarming
- Alarms should be visible at minimum of 2 m and easily audible (70 dB)
- All blood alarms should isolate the patient
 - Shut off the blood pump
 - Clamp the venous return
 - Stop UF



The blood side



Blood is pumped in the circuit @ 200-600 ml/min



The blood pump has 2 rollers compressing the tubing, forcing the blood along the tube
-adaptable to different sized tubing if required clinically



-can be operated manually in case of power loss

-calibrated to measure blood flow based on the internal diameter of the tubing

$BFR = RPM \times \text{tubing volume} (\pi \cdot r^2 \times \text{length})$


Possible sites for Air entry and blood loss






Excessively Negative Arterial pressure


- Low BP
- Catheter /needle malposition
- Kink/Clot
- Suction/spasm of vessel wall
- Stenosis in AVF
- Long small bore needles





The Venous Chamber

1. Allows for determining Positive pressure on the blood access site
2. Allows for separation of air bubbles (*excessive air, when present, may still enter the patient*)



High venous pressure alarm

- High Q_b , narrow needle
- Clot in the venous drip chamber, needle, venous limb
- Spasm of venous limb of AV access
- Kink in the venous return line
- Needle malposition

Low venous pressure alarm

Question

- 2 hours into HD, a patient c/o feeling weak. He has been c/o cramps for last 2 weeks.
- BP is 110/70 mm Hg (initial BP 150/95 mm Hg). HR 75/min, regular. Qb 500 ml/min. Qd 800 ml/min.
- Nurse checks HD machine, turns UF off, reduces blood flow to 200 ml/min
- HD is resumed



What happened? Venous Needle Dislodgement

Why did the venous alarm not sound?


Safety is the freedom from unacceptable risk
International Standard Organization 2007

The venous pressure reading during HD is determined by

- A. Intra access pressure
- B. Needle/Cannula geometry
- C. Blood flow rate and viscosity
- D. All of the above

Venous monitor: Pressures and Heights

- The fistula pressure
- Q_b , viscosity
- Flow resistance
- The height difference between AVF and the level in the venous chamber



Ribitsch W. et al Seminars in Dialysis 2013


In absence of blood flow

- Venous pressure reading = Intra access pressure
 - (when corrected for height difference between the access and the pressure monitor)
- ❖ When needle slips out from access but remains at the same height, the venous pressure will decrease by the amount of the access pressure

Needles may contribute significantly to venous pressure readings

□ Venous pressure


- Intra-access pressure
- Q_b
- Needle length
thickness
- Blood viscosity



With narrow gauge needles and high Q_b , the relatively higher flow resistance within needle may prevent the venous pressure to fall low enough to set off the alarm

The lower limit of the venous pressure alarms are usually set 30-40 mm below the access pressure

For the alarm to sound after VND, the pressure drop should exceed the pressure difference between the actual venous pressure and the lower alarm limit venous sensor



Fistulae have much lower pressure than grafts

AVG 60 mm Hg ; AVF 32 mm Hg

Over 90% AVG have access pressures >40

Only 30% AVF have access pressures >40




Fig. 2. Intra-access pressure (mmHg) in fistulae and grafts. VND: venous needle dislodgement

ates the lower detection limit of the

Lower detection limit
of venous sensor

Risk factors associated with venous needle dislodgement



- Never set wide alarm windows for VP monitors
- Check access and tubing when alarms are reset
- Secure tubing connections and cannula position
- Avoid covering access sites

Loss of 15-20% of blood volume over half an hour will cause death of a person with impaired sympathetic reflexes. **Guyton AC. 1991**

In just 5-7 minutes a patient on HD can lose **40%** of blood Volume from a **Venous Needle dislodgement**

UK Renal unit Survey [Clinical Directors and Lead renal nurses]


- Estimated prevalence/incidence of dislodgement
 - UK ~ 100/year (0-4 episodes/unit/year)
 - ~ **1/100,000 haemodialysis sessions**
- Severity
 - 1 death (0.6%)
 - 6.4% Moderate/Severe harm (e.g. hospitalization)
 - 93.0% No/Mild Harm

Many nurses and doctors are unaware of the limitations of venous alarms


Vascular Access Monitoring To Prevent Blood Leakage



Vascular access monitoring with blood sensor connected via Bluetooth to HD monitor



Gutter of drainage with blood sensor stopping blood pump






- Ultrasonic
- Located on return line
- Alarm activated by air in blood/saline or both


UABD lies below the bubble trap to prevent air passage

- Incidence of major air embolus 1:2000 treatments
- Usual volume required= 60-125 ml, esp. if injected rapidly

Air Microbubbles are found in Hemodialysis Patients




Microscopic finding of microbubbles of air (arrow) in pulmonary capillary; also, increased fibrosis of tissue (red areas).




Microscopic finding of microbubbles of air in the brain (arrow).

Microemboli in brain in HD




Microemboli can enter into HD patients without triggering alarms




- 54 HD patients (16 CVC, 38 AVF)
- Microemboli measured at AVF and Common Carotid artery before and during HD
- Significant increase in microembolic signals detected at both sites during HD

Preventing air embolism




Handling	Measure
Negative Pressure Connections	Connection between needle/CVC and arterial tube closed Heparin tube connector closed Arterial Pressure tube tightly applied
Priming of EC circuit	More volume is better (1.5l) Turn and knock on steam sterilized dialyzers
Blood Pump	Qb <200 ml/min
Venous air trap	Avoid low fluid level
Syringes	Take care to empty air
Dialyzer	Wet stored dialyzer preferred; careful connection to avoid air contamination
Air trap and dialysis device	Each brand has different conditions. Check with provider
EC components (tubing sets, dialyzers)	Each brand has different conditions. Check with provider

The Dialysate pathway



Simplified Fluid Pathway




Adapted from Handbook of Dialysis Therapy. Nissenson and Fine. 4th edition,
Saunders Elsevier Pubs.

The three stream method of preparing a bicarbonate based dialysate

Organic acid

Glacial acetic acid/sodium diacetate/
citric acid



45 y/ white male on HD X 4 y.

Cause of ESRD: Chronic GN

Pre Labs (midweek): K 5, Na 138, CO₂ 30, Cl 112, Ca 9.3, Phosphorus 6

Dialysis Prescription

Pre 86.5 kg

Target 84 kg

Temp 37.2 c

Dialysate Temp 36 c

Qb 350

Qd 800

Time 240 min

Dialysate 3K/2.5 Ca/ 39

Bicarb


Pre K 5, Na 138, CO₂

30, Cl 112, Ca 9.3,

Phosphorus 6

Dialyzer F80

2 hours into dialysis, patient c/o palpitations.
HR 150, irr.irr, BP 90 mm systolic, RR 24/min,
SaO₂ 85% on room air



- ❖ Patient is disconnected and HD is stopped.
- ❖ Oxygen, IV Saline
- ❖ Stat blood draw: K 2.0 meq/L, HCO₃ 36 meq/L
- ❖ IV Potassium, Cardiology consult, CICU
- ❖ Dialysate and Patient Chemistry match analysis

The final dialysate solution

Na Bicarbonate (37 mEq/L) + Acetic acid (4 Mm)

Carbonic Acid + Na acetate

Carbon dioxide + water + Na acetate (4 mEq/L)

The loss of bicarbonate is balanced by a gain in acetate
(bicarbonate precursor)

Total buffer base = Na Bicarbonate (33 mEq/L) + Na acetate (4 mEq/L)
= 37 mEq/L

It is vital to know the total buffer base available for administration to the patient via the final dialysate at a specific machine setting and concentrate pairing

The dialysate acid concentrate can contain acetic acid, acetate or citrate, in amounts ranging from 1.5 to 8 mEq/L

This generates bicarbonate in the body!

- Possible in all currently marketed dialysate concentrate products containing acetate, acetic acid, or citrate.
- 50 HD patients hospitalized in October 2010
- Their outpatient dialysate prescription included a **39 mEq/L bicarbonate solution and an acid concentrate which contained 8 mEq/L of acetate (total bicarbonate of 43 mEq/L)**.
- At presentation, the patients' mean serum bicarbonate level was **31.3 mEq/L and 54 percent had a serum bicarbonate >30 mEq/L**.

Pande S, Raja R, Bloom E et al. Effect of dialysate baths on serum bicarbonate levels in hemodialysis patients. AJKD 2011; 57(4): A75 (Abstract #234))



Medical Devices



Home | Medical Devices | Medical Device Safety | Medical Device Recalls

Sodium Diacetate, has equal parts of acetic acid and sodium acetate.

The product dialysate has a total of 8 mEq/L of acetate, 4 from each component.


With a starting Bicarbonate concentration of 37 mEq/L, the total buffer base in the final dialysate will be 41 mEq/L
 $(37 - 4 = 33 + 4 + 4 = 41)$

(acetic acid plus acetate) contained in Fresenius' Naturalyte Liquid and Granule Dry Acid Concentrate. Inappropriate prescription of these products can lead to a high serum bicarbonate level in patients undergoing hemodialysis. This may contribute to metabolic alkalosis, which is a significant risk factor associated with low blood pressure, hypokalemia, hypoxemia, hypercapnia and cardiac arrhythmia, which, if not appropriately treated, may culminate in cardiopulmonary arrest. This product may cause serious adverse health consequences, including death.

FDA has issued a general safety communication related to inappropriate prescription and resultant alkali dosing errors in the dialysate concentrates used in hemodialysis.


Proper mixing monitored by conductivity monitors [normal range 12-16 mS/cm]

- Determined by ionic constituents
- Greater the number of ions greater the conductivity
- Should be routinely checked by measuring dialysate Na concentration




Misra M. The basics of HD equipment.
Hemodialysis International 2005; 9: 30-36

Low conductivity Alarm



High conductivity Alarm




Adapted from Handbook of Dialysis Therapy. Nissenson and Fine. 4th edition,
Saunders Elsevier Pubs.

Potential proportioning problems

- Causes
 - Wrong concentrate (note color coding: **Red** for acid, **Blue** for base)
 - Poor mixing
 - Clogged filters
 - Crystallization in the system
 - Human disarming of the switches
- Outcome
 - ❖ High or low plasma sodium
 - ❖ High or low plasma potassium
 - ❖ High calcium/magnesium
 - ❖ High or low plasma osmolality

Simplified Fluid Pathway



Adapted from Handbook of Dialysis Therapy. Nissenson and Fine. 4th edition,
Saunders Elsevier Pubs.

Dialysate Alarms

Alarms should interrupt the supply of dialysate

- Check for ‘no flow’ and lack of dialysate stream at the dialyzer

Problem Solving

- Is the concentrate container empty? **[low conductivity]**
- Is the water inlet pressure normal? **[high conductivity]**
- Are there any water leaks/puddles beneath the mixing chambers? **[high conductivity]**
- Never adjust conductivity setting when the patient is on dialysis

Conclusions

- Hemodialysis machine is part of an integrated system.
- Safety alarms and monitors are not foolproof
- User errors are common
- Basic understanding of functions of blood and dialysate side circuits and equipment helps avoid complications
- Detailed technical knowledge is not necessary but fundamental knowledge must be acquired
- Patient safety is paramount!